

Bath Step - User Instructions


(Model: AB0028)

MI0247 issue 1

The Bath Step by Mangar International is strong and lightweight and the textured surface ensures maximum security when stepping into and out of a bath. The height adjusting feet provide a step height of 15cm or 10cm.

Contents:

- 1 x Bath Step platform (A) with detachable feet (B).
- 2 x Height adjusting feet (C).




How to use the Bath Step

- Set the height to either 15cm or 10cm by fitting or removing the height adjusting feet (C).

WARNING: DO NOT use the step without the detachable feet (B) fitted.

- Ensure the Bath Step is dry and free from soap or oil residue.
- Position the Bath Step close to the bath side.
- Ensure that the Bath Step is stable and stand on the step with both feet centrally on the platform.
- Step directly into the bath.




Cleaning

Clean the Bath Step with a non-scratch detergent.

Specifications:

Maximum user weight:
200Kg (31st)



Warranty

Mangar International guarantees this product for a period of one year.

The guarantee is given against defects in the product. The guarantee is not valid in the case of fair wear and tear, wilful damage, common neglect, misuse or unauthorised alteration or repair by unqualified persons.

The following statements do not affect the consumer's statutory rights. This product is sold on the understanding that in the event of any defect of manufacture or material appearing within one year of the date of receipt of this product, the product will be repaired or replaced free of charge providing that:

- Reasonable evidence is provided (e.g. purchase invoice, delivery note) that the product was not purchased more than one year prior to the date of the claim.
- The defect is not attributable to accidental damage (either in transit or otherwise), misuse or unauthorised repair.

Mangar International reserves the right to alter product specifications and/or any of the information contained within this document without notice.



Mangar International

Presteigne, Powys, Wales, UK, LD8 2UF

Tel: +44 (0)1544 267674

Fax: +44 (0)1544 260287

customercare@mangar.co.uk

www.mangarinternational.co.uk

