



Max: 150 kg/330 lbs

OneWaySlide is used for repositioning farther back in a chair/wheelchair. OneWaySlide is designed to glide easily backward, but not forward. It is placed so that the patient remains sitting on it.

OneWaySlide is part of a series of different models of transfer assistive devices. The choice of model and material is determined by the patient's functional capacity and the situations in which the product is to be used.

RoMedic's products together form the SystemRoMedic concept, which takes a holistic approach to patient transfers. SystemRoMedic has four categories of transfers: transfer, positioning, support and lifting.



Functional inspection

Visual inspection

Check the condition and function of the product regularly, and always after laundering. Check to ensure that the fabric is not worn, damaged or bleached. If there are signs of wear, the product must be discarded.




Always read the instructions

Read instructions for all assistive devices used in connection with transfers.

To prevent discomfort, exercise care and caution when using the product.



Functional inspection



Test OneWaySlide's friction by rolling it and noting how easily it glides. Check to ensure that it resists movement in the opposite direction. If possible, compare with a new OneWaySlide.

Placement of OneWaySlide



Check to ensure that the OneWaySlide lies securely on the seat cushion. A rough material on the top of the seat cushion gives the best result.

Place the OneWaySlide as indicated by the arrow on the product label and allow it to hang about 10 cm over the front edge of the cushion.

NOTE! Place OneWaySlide according to instructions on product label. Incorrect placement may entail a risk of the patient sliding out of the chair.



Using OneWaySlide



If the patient slides forward in the wheelchair, the caregiver can help them to lean forward and provide counter pressure, thereby enabling the patient to slide back in the chair.

Make sure the knees are raised to the same level as the hips. This can be done with the aid of a FootStool.

Material:

Polyester, polyamide, cotton

Size

2069 50 x 40

2070 50 x 45

Care of the product:

Read the product label.



Do not use rinsing agent. To ensure maximum material life, hang dry or tumble dry on low heat.

When in doubt, contact your RoMedic representative for product trials, advice or information.

If you require additional information, contact your transfer specialist or our representative, or visit www.romedic.com. Make sure you always have the most recent version of the manual. Manuals are available for downloading from our website;

<http://www.romedic.com>

RoMedic AB

Box 640

SE-175 27 Järfälla

SWEDEN

Tel: +46 (0)8-594 773 33

Fax: +46 (0)8-594 773 44

Internet: www.romedic.com

E-mail: office@romedic.com

RoMedic[®]
A company in the Handicare Group