


Code	Overall Length	Length between Centres	Handle Length	Finish
GR12C	14½" (37cm)	12" (31cm)	7" (18cm)	White
GR12W	14¾" (37cm)	12" (31cm)	7" (18cm)	Chrome
GR16C	18½" (47cm)	16" (41cm)	11" (28cm)	White
GR16W	18½" (47cm)	16" (41cm)	11" (28cm)	Chrome
GR18C	20¾" (51cm)	18" (46cm)	14" (36cm)	White
GR18W	20¾" (51cm)	18" (46cm)	14" (36cm)	Chrome
GR32C	34½" (88cm)	32" (81cm)	27" (69cm)	White
GR32W	34½" (88cm)	32" (81cm)	27" (69cm)	Chrome

DIAGRAMS ARE SUPPLIED FOR INFORMATION ONLY.  
THE MODEL RECEIVED MAY DIFFER SLIGHTLY FROM THE IMAGE.



Ainley's Industrial Estate, Eiland HX5 9JP  
T: +44 (0) 1422 314488  
E: enquiries@drivemedical.co.uk  
W: www.drivemedical.co.uk/technical

NOTES:

TITLE:

## GRAB RAIL DIMENSIONS

DATE: 08/05/13

SHEET No: 1 of 1

DRAWN: PMK

CHECKED:

REVISION No: —

SCALE: NTS

DWG No: GR