

Why choose any other FIXED bathlift?



AquaJoy Bathlifts Ltd
Consett Business Park
Consett
County Durham
DH8 6BN
United Kingdom

Tel: +44 (0)1207 501555
Fax: +44 (0)1207 599789
email: info@aqua-joy.com
www.aqua-joy.com



NEW

SAVE space

SAVE money

SAVE time and effort

SAVE back strain

OUT NOW!

AquaJoy Saver & Junior Saver Bathlifts

So why choose any other **FIXED** bathlift?

SUITABLE FOR ADULTS AND CHILDREN

- Adaptable accessories for users with varying ages and needs
- BARIATRIC WEIGHT LIMIT AS STANDARD 160KG (24 STONE)**

- Suitable for heavier clients

FITS ALMOST ANY STYLE OF BATH

- Ideal for deep baths – maximum seat height 465mm
- Ideal for shallow baths – minimum seat height 78mm

MAXIMISES LEGROOM IN THE BATH

- Bathlift seat locates as far back as possible in bath
- Short footprint and seat

INTERCHANGEABLE SPARE PARTS AND ACCESSORIES

- Saving time, effort, stockholding and money!



Moving and Handling – made easy!

HEAVIEST PART ONLY 4.65KGS

- Reduces the risk of moving and handling injuries
- Reduces back strain when fitting and removing

COMPACT DESIGN WHICH CURVES AROUND THE BODY

- Follows guidelines for the safe lifting of loads
- Allows an unobstructed view

EASY RELEASE SUCTION FEET WITH LIFTING BAR AS STANDARD

- Limits stooping and stretching when lifting the Bathlift in and out of the bath
- No struggling under the bathlift seat to detach individual suction feet



Easy to clean and maintain

CAN BE USED WITH OR WITHOUT COVERMATS

- Detachable and can be washed in a standard washing machine up to 60 degrees protecting against unwanted germs and bacteria

SMOOTH, CREVICE FREE SURFACE

- Dirt and used water do not get trapped
- Makes cleaning quick and easy

STREAMLINED FRAMEWORK

- Provides open access to the Bathlift and bath for cleaning when in situ

Simple and safe to operate

LIGHTWEIGHT, FULLY WATERPROOF, FLOATING HANDCONTROL

- Comfortable to hold and operate
- Suction pads on reverse to stick to the side of the bath or to tiles for ease of use

CONVEX AND CONCAVE BUTTONS

- Soft to touch
- Easy to operate for users with limited dexterity or visual impairment

SIMPLE CHARGING PROCESS

- Pronged connectors for ease of use

BATTERY INDICATOR LIGHT

- Indicates when recharging is required

CHARGER INDICATOR LIGHT

- Indicates charging is in progress

SAFETY FEATURE

- The Bathlift will not lower unless there is sufficient charge to raise the client back to the top of the bath

