

# POWERSTROLL

Wheelchair Accessories

**Power**STROLL

S•Drive

*“Takes the strain  
out of pushing  
a wheelchair!”*



CE

This excellent Powerstroll will turn most manual wheelchairs into carer controlled powered wheelchairs. It takes the strain out of pushing a heavy wheelchair and makes easy work of hills, ramps and uneven surfaces. Its unique one click docking station (patented) means that it can be fitted or removed easily in seconds.

**drive**<sup>TM</sup>  
[www.drivemedical.co.uk](http://www.drivemedical.co.uk)

The Powerstroll can easily be fitted and removed within seconds through its unique one-click docking station and the new push and slide battery removal system.



Easy to fit & remove

## Superb Comfort, Style & Performance...

- » Allows speed of up to 4mph\*
  - » Travels up to 10 miles (16 km) on a full battery charge\*
  - » The Powerstroll can easily be fitted or removed within seconds
  - » The excellent P G Drives S-Drive controller with fully programmability & comprehensive fault management system
  - » New push and slide battery removal system
  - » It is very portable, making transportation simple
  - » Fits most standard wheelchairs up to 20" seat width
  - » Additional parts available for wider wheelchairs
  - » Hand control with variable speed dial, forward/reverse function, illuminated battery gauge and fingertip trigger control
  - » Unique one-click docking station for ease of assembly (patented)
  - » Comes complete with a Carry Holdall with shoulder strap
  - » Simple to manoeuvre and reverse
  - » The Powerstroll has a powerful 150 watt (200 watt-HD), 24V motor, capable of coping with most terrains
  - » Maintenance free airline friendly 24V 12ah valve regulated sealed lead acid battery pack
  - » 2A charger as standard for off-board charging
  - » Weight capacity 135kg (21st), this includes the wheelchair and user combined
  - » Heavy Duty model - weight capacity 180kg (28st)
- \* Subject to conditions of use



Hand control with variable speed dial, forward/reverse function, illuminated battery gauge and fingertip trigger control.

CODE: **PWCPP009 & PWCPP009HD**  
DESCRIPTION: **Powerstroll S Drive**



Total Weight:  
19.2kg (42lb)



Weight Without Battery:  
10.2kg (22.5lb)



Weight Battery Pack Only:  
9kg (19.5lb)



Max Range:  
10 Miles (16km)\*



Max Gradient:  
8°



Max Weight Capacity:  
135kg (21st) Standard  
180kg (28st) HD Model



Fits most standard wheelchairs up to 20" seat width

Leaflet Code: **LL144**

If you would like to view this product please contact your local stockist, details below.